



Product Specification

Sample changer ATD type GS 301, for gas sampler GS 301 (item no 90.170.350) & tubes 6 x 89 mm

Page 1



Product description

Order number 92.170.354
Product description Sample changer ATD type GS 301, for gas sampler GS 301 (item no 90.170.350) & tubes 6 x 89 mm, 1 piece(s)/case

Product characteristics

Suitable for gas sampler GS 301 (item no 90.170.350) & tubes 6 x 89 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 1
Case volume 0.0011 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-12-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Sample changer ATD type GS 301, for gas
sampler GS 301 (item no 90.170.350) & tubes 6 x
89 mm

Page 2

Weight of product	0.55 kg
Weight of case	0.7 kg
EAN of inner box	4038917101219
EAN case	4038917051538

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-12-04

This document was prepared by EDP
support and is valid without signature.



SARSTEDT