



Product description

Order number 92.170.355
Product description Sample changer GE type GS 301, for gas sampler GS 301 (item no 90.170.350) & tubes 6 x 178 mm, 1 piece(s)/case

Product characteristics

Suitable for gas sampler GS 301 (item no 90.170.350) & tubes 6 x 178 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging case
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 1

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification

**Sample changer GE type GS 301, for gas sampler
GS 301 (item no 90.170.350) & tubes 6 x 178 mm**

Page 2

Case volume	0.0016 cbm
Weight of product	0.6 kg
Weight of case	0.8 kg
EAN of inner box	4038917101226
EAN case	4038917051545

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-05-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT