



## Product description

Order number	92.2620.304
Product description	length: 1.5 m, material: silicone, transparent, for IVARO, 1 piece(s)/bag

## Product characteristics

Suitable for	IVARO
--------------	-------

## Size

Diameter	0.8 mm
Length of product	1.5 m

## Material & colours

Product material	silicone
Colour of product	transparent

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	10
Depth of case	280 mm
Width of case	220 mm
Height of case	100 mm
Case volume	0.0062 cbm
Weight of product	0.0005 kg
EAN of inner box	4038917584111
EAN case	4038917584104

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

