



Product Specification
Gas sampler GS 212, incl. Mains adapter GN
100?240/15 V

Page 1



Product description

Order number	90.170.310
Product description	Gas sampler, GS 212, incl. Mains adapter GN 100?240/15 V, feed rate: 0.2 - 12 l/min, feed volume: 1 - 9,999 l

Product characteristics

Power consumption during operation	24 W
---------------------------------------	------

Size

Width of product	410 mm
Length of product	220 mm

Purity & certification

Product category	no medical device IVD
CE certified	CE
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Gas sampler GS 212, incl. Mains adapter GN
100?240/15 V

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.06 cbm
Weight of product	11 kg
Weight of case	13.24 kg
EAN of inner box	4038917101141
EAN case	4038917051484

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

