



### **Product description**

Order number	92.170.390
Product description	Transport case GS 100, for Microbial air samplers, 1 piece(s)/case

### **Product characteristics**

Suitable for	Microbial air samplers
--------------	------------------------

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	490 mm
Width of box	420 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Height of box	295 mm
Case volume	0.0286 cbm
Weight of product	2.7395 kg
Weight of case	2.74 kg
EAN of inner box	4038917101295
EAN case	4038917051590

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

