



Product description

Order number	92.170.390
Product description	Transport case GS 100, for Microbial air samplers, 1 piece(s)/case

Product characteristics

Suitable for	Microbial air samplers
--------------	------------------------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0286 cbm
Weight of product	2,7395 kg
Weight of case	2,74 kg
EAN of inner box	4038917101295

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN case

4038917051590

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT