



### Product description

Order number	80.626
Product description	Forensic swab, round, ETO treated, ISO 18385, stem length: 150 mm, stem material: wood, swab tip round, swab material: cotton, 1 piece(s)/blister

### Product characteristics

Application	Removal of forensic samples
-------------	-----------------------------

#### Size

Diameter	2.2 mm
Width of product	2.2 mm
Length of product	150 mm
Length of stem	150 mm

### Material & colours

Colour of product	multi-coloured
Swab material	cotton
Stem material	wood

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Satisfies the requirement	ISO 18385
Product category	not regulated
CE certified	no CE mark
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	68000
Depth of box	200 mm
Width of box	105 mm
Height of box	75 mm
Depth of case	428 mm
Width of case	398 mm
Height of case	135 mm
Case volume	0.023 cbm
Weight of product	0.0012 kg
Weight of case	1.95 kg
EAN of inner box	4038917279345
EAN case	4038917279338

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

