



### Product description

|                     |   |
|---------------------|---|
| Order number        | 80.626  |
| Product description | Forensic swab, round, ETO treated, ISO 18385, stem length: 150 mm, stem material: wood, swab tip round, swab material: cotton, 1 piece(s)/blister |

### Product characteristics

|             |                             |
|-------------|-----------------------------|
| Application | Removal of forensic samples |
|-------------|-----------------------------|

#### Size

|                   |        |
|-------------------|--------|
| Diameter          | 2.2 mm |
| Width of product  | 2.2 mm |
| Length of product | 150 mm |
| Length of stem    | 150 mm |

### Material & colours

|                   |                |
|-------------------|----------------|
| Colour of product | multi-coloured |
| Swab material     | cotton         |
| Stem material     | wood           |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | ISO 18385               |
| Product category          | no medical device   IVD |
| Purity standard           | ETO treated             |
| Sterilisation             | Gas sterilisation       |
| Batched                   | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 68000         |
| Depth of case                 | 428 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 135 mm        |
| Case volume                   | 0.023 cbm     |
| Weight of product             | 0.0012 kg     |
| Weight of case                | 1.95 kg       |
| EAN of inner box              | 4038917279345 |
| EAN case                      | 4038917279338 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

