



### **Product description**

Order number	80.626
Product description	Forensic swab, ETO treated, ISO 18385, passive drying, stem length: 150 mm, stem material: wood, swab material: cotton, 1 piece(s)/blister

### **Product characteristics**

Application	Removal of forensic samples
-------------	-----------------------------

#### **Size**

Diameter	2,2 mm
Width of product	2,2 mm
Length of product	150 mm
Length of stem	150 mm

#### **Material & colours**

Colour of product	multi-coloured
Swab material	cotton
Stem material	wood

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Satisfies the requirement	ISO 18385
Product in accordance with MDD or IVDD	no medical device   IVD
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Lot or serial number	yes

### **Packaging**

Minimum order qty.	1000
Type of inner packaging	blister
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	68000
Depth of case	428 mm
Width of case	398 mm
Height of case	135 mm
Case volume	0,023 cbm
Weight of product	0,0012 kg
Weight of case	1,95 kg
EAN of inner box	4038917279345
EAN case	4038917279338

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

