



Product description

Order number	80.630
Product description	Forensic swab, round, in the tube with ventilation membrane, ETO treated, ISO 18385, passive drying, stem length: 85 mm, stem material: PS, swab tip round, swab material: viscose, tube: PP, 1 piece(s)/transport tube

Product characteristics

Label/ Print	with paper label
Cap	standard cap
Type of cap	slip-on lid
Design	in the tube with ventilation membrane
Application	Removal of forensic samples

Size

Length of product	175 mm
Length excluding cap	105 mm
Length of stem	85 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Colour of product	white
Cap material	Polypropylene (PP)
Colour of cap	white
Swab material	viscose
Stem material	Polystyrene (PS)
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	ISO 18385
Product category	no medical device IVD
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	transport tube
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	33000
Depth of box	181 mm
Width of box	97 mm
Height of box	98 mm
Depth of case	498 mm
Width of case	373 mm
Height of case	115 mm
Case volume	0,0214 cbm
Weight of product	0,0074 kg
Weight of case	4,52 kg
EAN of inner box	4038917275026
EAN case	4038917275019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

