



Product description

Order number	82.1969.002
Product description	Deep Well MegaBlock® , 96 well, 0.5 ml, PP, round base, free from DNA, DNase/RNase, pyrogens/endotoxins, material: PP, SLAS standard, 7 piece(s)/bag

Product characteristics

Type of cap	film
Design	SLAS standard
Application	collection and storage of diagnostic specimens
Base shape	round base
Number of wells	96
Volume of work	1 µl - 0,5 ml
Centrifugation max (RCF)	2000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	0,5 ml
Diameter	7,2 mm
Width of product	85,1 mm
Height of product	28,5 mm
Length of product	127,3 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNA, DNase/RNase, pyrogens/endotoxins
Batched	yes

Packaging

Minimum order qty.	56
Type of inner packaging	bag
Piece(s) / inner box	7
Piece(s) / case	56
Piece(s) / pallet	4480
Case volume	0,0181 cbm
Weight of product	0,057 kg
Weight of case	3,74 kg
EAN of inner box	4038917346320
EAN case	4038917346313

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

