



Product description

Order number	82.1972.800
Product description	Tip comb, compatible with KingFisher? Flex/Presto, Bio Sprint 96, free from DNA/DNase/RNase/PCR inhibitors, material: PP, 2 piece(s)/bag

Product characteristics

Suitable for	KingFisher? Flex/Presto, Bio Sprint 96
Autoclavable (max. °C)	121 °C

Size

Diameter	7.52 mm
Width of product	85 mm
Height of product	41.85 mm
Length of product	127.76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tip comb, free from DNA/DNase/RNase/PCR inhibitors, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase/PCR inhibitors
Batched	yes

Packaging

Minimum order qty.	20
Type of smallest subpackaging	bag
Piece(s) / inner box	2
Piece (s) / outer case	20
Piece(s) / pallet	2160
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0543 kg
Weight of case	1.28 kg
EAN of inner box	4038917958080
EAN case	4038917958097

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-12

This document was prepared by EDP support and is valid without signature.



SARSTEDT