



Product description

Order number	95.2093
Product description	Film, free from DNase/RNase, transparent, pierceable, application area: -80 - 120 °C

Product characteristics

Type of cap	film
Light protection	no
Suitable for qPCR	no
Minimum temperature	-80 °C
Maximum temperature	120 °C

Size

Width of product	75 mm
Length of product	118 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Film, free from DNase/RNase, transparent,
pierceable**

Page 2

Material & colours

Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard free from DNase/RNase
Batched yes

Packaging

Minimum order qty. 20
Piece(s) / inner box 20
Piece(s) / case 240
Case volume 0,0142 cbm
Weight of product 0,012 kg
Weight of case 3 kg
EAN of inner box 4038917335706
EAN case 4038917335690

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT