



**Product Specification**  
**Disposal bags, 2 l, (LxW): 300 x 200 mm, PP,**  
**transparent, without print**

**Page 1**



### Product description

Order number	86.1197
Product description	Disposal bags, volume: 2 l, (LxW): 300 x 200 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, without print, 100 piece(s)/dispenser

### Product characteristics

Label/ Print	without print
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

### Size

Width of product	200 mm
Length of product	300 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-10-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Safety data sheet	no
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	dispenser
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	63000
Depth of box	73 mm
Width of box	73 mm
Height of box	211 mm
Depth of case	388 mm
Width of case	253 mm
Height of case	165 mm
Case volume	0.0162 cbm
Weight of product	0.0056 kg
Weight of case	6.1 kg
EAN of inner box	4038917085717
EAN case	4038917004343

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

