



**Product description**

Order number	86.1198
Product description	Disposal bags, volume: 7 l, (LxW): 500 x 300 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, without print, 50 piece(s)/bag

**Product characteristics**

Label/ Print	without print
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

**Size**

Width of product	300 mm
Length of product	500 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Disposal bags, 7 l, (LxW): 500 x 300 mm, PP,**  
**transparent, without print**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Safety data sheet	no
Batched	no

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	48000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0135 kg
Weight of case	7 kg
EAN of inner box	4038917085724
EAN case	4038917004350

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

