



### **Product description**

Order number	86.1199
Product description	Disposal bags, volume: 24 l, (LxW): 780 x 400 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, without print, 50 piece(s)/bag

### **Product characteristics**

Label/ Print	without print
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

### **Size**

Width of product	400 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal bags, 24 l, (LxW): 780 x 400 mm, PP,**  
**transparent, without print**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Latex-free yes  
Safety data sheet no  
Batched no

**Packaging**

Minimum order qty. 250  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 250  
Piece(s) / pallet 16500  
Depth of case 438 mm  
Width of case 268 mm  
Height of case 155 mm  
Case volume 0.0182 cbm  
Weight of product 0.0271 kg  
Weight of case 7.05 kg  
EAN of inner box 4038917085731  
EAN case 4038917004367

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-01-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**