



Product description

Order number	86.1201
Product description	Disposal bags, volume: 7 l, (LxW): 500 x 300 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	red
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

Size

Width of product	300 mm
Length of product	500 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal bags, 7 l, (LxW): 500 x 300 mm, PP,
transparent, with print Biohazard

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Safety data sheet no
Batched no

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 25000
Depth of box 305 mm
Width of box 135 mm
Height of box 25 mm
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.0141 kg
Weight of case 7.3 kg
EAN of inner box 4038917085755
EAN case 4038917004381

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-03-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT