



Product Specification
Disposal bags, 7 l, (LxW): 500 x 300 mm, PP,
transparent, with print Biohazard



Product description

Order number	86.1201
Product description	Disposal bags, volume: 7 l, (LxW): 500 x 300 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	red
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

Size

Width of product	300 mm
Length of product	500 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal bags, 7 l, (LxW): 500 x 300 mm, PP,
transparent, with print Biohazard

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	25000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0141 kg
Weight of case	7.3 kg
EAN of inner box	4038917085755
EAN case	4038917004381

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

