



Product Specification
Disposal bags, 24 l, (LxW): 780 x 400 mm, PP,
yellow, with print Biohazard

Page 1



Product description

Order number	86.1202.103
Product description	Disposal bags, volume: 24 l, (LxW): 780 x 400 mm, material: PP, yellow, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: blue, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	blue
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

Size

Width of product	400 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal bags, 24 l, (LxW): 780 x 400 mm, PP,
yellow, with print Biohazard

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	9000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0289 kg
Weight of case	7.5 kg
EAN of inner box	4038917160438
EAN case	4038917160421

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT