



**Product Specification**  
**Disposal bags, (LxW): 900 x 320 mm, LD-PE,**  
**yellow, Biohazard**

**Page 1**



### Product description

Order number	86.9925.429
Product description	Disposal bags, (LxW): 900 x 320 mm, material: LD-PE, yellow, film thickness: 50 µm, print Biohazard, 50 piece(s)/bag

### Size

Width of product	320 mm
Length of product	900 mm

### Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	yellow

### Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Disposal bags, (LxW): 900 x 320 mm, LD-PE,**  
**yellow, Biohazard**

**Page 2**

**Packaging**

Minimum order qty.	150
Type of inner packaging	bag
Piece(s) / inner box	50
Piece(s) / case	150
Piece(s) / pallet	10000
Depth of case	428 mm
Width of case	358 mm
Height of case	234 mm
Case volume	0.0357 cbm
Weight of case	13.62 kg
EAN of inner box	4038917155564
EAN case	4038917400374

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**