



Product description

Order number	95.1297
Product description	Table stand, for destruction bags 300 x 200 mm (item no 86.1197), material: epoxy resin, white, (HxØ level1/Ø level2): 250 x 155/125 mm, incl. dispensing package with 100 bags, 1 piece(s)/case

Product characteristics

Suitable for	destruction bags 300 x 200 mm (item no 86.1197)
Accessories	geeignet für Vernichtungsbeutel 200 x 300 mm (Artikel-Nr. 86.1197)

Size

Height of product	250 mm
-------------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification

Table stand, for destruction bags 300 x 200 mm
(item no 86.1197), epoxy resin, white

Page 2

Material & colours

Product material epoxy resin
Colour of product white

Purity & certification

Product category no medical device | IVD
Latex-free yes
Safety data sheet no
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 108
Depth of box 73 mm
Width of box 73 mm
Height of box 211 mm
Depth of case 308 mm
Width of case 223 mm
Height of case 195 mm
Case volume 0.0134 cbm
Weight of product 0.658 kg
Weight of case 0.88 kg
EAN of inner box 4038917086073
EAN case 4038917004664

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT