



Product Specification
Table stand, for destruction bags 300 x 200 mm
(item no 86.1197), epoxy resin, white

Page 1



Product description

Order number	95.1297
Product description	Table stand, for destruction bags 300 x 200 mm (item no 86.1197), material: epoxy resin, white, (HxØ level1/Ø level2): 250 x 155/125 mm, incl. dispensing package with 100 bags, 1 piece(s)/case

Product characteristics

Suitable for	destruction bags 300 x 200 mm (item no 86.1197)
Accessories	geeignet für Vernichtungsbeutel 200 x 300 mm (Artikel-Nr. 86.1197)

Size

Height of product	250 mm
-------------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Table stand, for destruction bags 300 x 200 mm
(item no 86.1197), epoxy resin, white

Material & colours

Product material	epoxy resin
Colour of product	white

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	108
Depth of box	73 mm
Width of box	73 mm
Height of box	211 mm
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.658 kg
Weight of case	0.88 kg
EAN of inner box	4038917086073
EAN case	4038917004664

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

