



### Product description

Order number	82.1581.120
Product description	ELISA plate, 96 well, flat base, without lid, material: PS, black, Medium Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

### Product characteristics

Label/ Print	with print
Base shape	flat base
Number of wells	96
Vessel type	Plate

### Size

Max. volume	390 µl
Width of product	85.2 mm
Height of product	14.35 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	black

**Purity & certification**

Product category	no medical device   IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	50
Piece(s) / pallet	2800
Depth of case	279 mm
Width of case	232 mm
Height of case	365 mm
Case volume	0.0236 cbm
Weight of product	0.0408 kg
Weight of case	2.86 kg
EAN of inner box	4038917365338
EAN case	4038917365321

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

