



### **Product description**

Order number	82.1582.001
Product description	Micro test plate, 96 well, slip-on lid, round base, lid assembled, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 1 piece(s)/blister

### **Product characteristics**

Label/ Print	with plastic label
Closure type	slip-on lid
Base shape	round base
Number of wells	96
Vessel type	Plate

### **Size**

Max. volume	310 µl
Width of product	85.2 mm
Height of product	16.55 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                    transparent

**Purity & certification**

Product category                    no medical device | IVD  
Purity standard                    free from DNA/DNase/RNase, free from  
   pyrogens/endotoxins, non-cytotoxic, sterile  
  
Sterilisation                        Gamma irradiation  
Batched                                yes

**Packaging**

Minimum order qty.                50  
Type of smallest                    blister  
subpackaging  
Piece(s) / inner box                1  
Piece (s) / outer case              50  
Piece(s) / pallet                    3200  
Depth of case                      332 mm  
Width of case                       297 mm  
Height of case                      215 mm  
Case volume                        0.0212 cbm  
Weight of product                  0.0498 kg  
Weight of case                      3.19 kg  
EAN of inner box                   4038917131063  
EAN case                            4038917347969

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

