



Product Specification
Micro test plate, 96 well, round base, PS,
transparent

Page 1



Product description

Order number	82.1582
Product description	Micro test plate, 96 well, round base, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Base shape	round base
Number of wells	96
Vessel type	Plate

Size

Max. volume	310 µl
Width of product	85.2 mm
Height of product	14.35 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	100
Piece(s) / pallet	3200
Depth of case	359 mm
Width of case	275 mm
Height of case	365 mm
Case volume	0.036 cbm
Weight of product	0.0343 kg
Weight of case	4.33 kg
EAN of inner box	4038917087100
EAN case	4038917355469

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.