



Product description

Order number	82.1583.001
Product description	Micro test plate, 96 well, slip-on lid, base shape: conical, lid assembled, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with plastic label
Closure type	slip-on lid
Base shape	conical
Number of wells	96
Vessel type	Plate

Size

Max. volume	290 µl
Width of product	85.2 mm
Height of product	16.55 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	3200
Depth of case	332 mm
Width of case	297 mm
Height of case	215 mm
Case volume	0.0212 cbm
Weight of product	0.0518 kg
Weight of case	3.29 kg
EAN of inner box	4038917131094
EAN case	4038917347976

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

