



### Product description

Order number	82.1584
Product description	Lid, for ELISA and microtest plates in 96-well format, material: PS, with condensation rings, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

### Product characteristics

Design	with condensation rings
Suitable for	ELISA and microtest plates in 96-well format

### Material & colours

Product material	Polystyrene (PS)
------------------	------------------

### Purity & certification

Product category	no medical device   IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Lid, for ELISA and microtest plates in 96-well  
format, PS, with condensation rings**

**Page 2**

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	100
Piece(s) / pallet	11200
Depth of case	282 mm
Width of case	187 mm
Height of case	235 mm
Case volume	0.0124 cbm
Weight of product	0.0153 kg
Weight of case	1.86 kg
EAN of inner box	4038917087124
EAN case	4038917005593

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-19

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**