




Product Specification
Micro test plate, Terasaki, 60 well, slip-on lid,
base shape: conical, PS, transparent

Page 1



Product description

Order number	83.9923.972
Product description	Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, lid assembled, material: PS, transparent, sterile, 10 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Base shape	conical
Number of wells	60
Vessel type	Plate
Type of surface treatment	TC-treated

Size

Max. volume	10 µl
Width of product	57.8 mm
Height of product	9.7 mm
Length of product	85.1 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Micro test plate, Terasaki, 60 well, slip-on lid,
base shape: conical, PS, transparent

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	650
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	650
Piece(s) / pallet	16250
Depth of case	608 mm
Width of case	324 mm
Height of case	264 mm
Case volume	0.052 cbm
Weight of product	0.0173 kg
Weight of case	12.2 kg
EAN of inner box	4038917154710
EAN case	4038917354363

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT