



Product Specification
**Weigh tray, 250 ml, (LxW): 128 x 128 mm, PVC,
white**

Page 1



Product description

Order number 71.9923.210
Product description Weigh tray, 250 ml, (LxW): 128 x 128 mm, material: PVC,
white, 100 piece(s)/bag

Size

Sample volume 250 ml
Max. volume 300 ml
Width of product 128 mm
Length of product 128 mm

Material & colours

Product material Polyvinylchloride (PVC)
Colour of product white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Weigh tray, 250 ml, (LxW): 128 x 128 mm, PVC,
white**

Page 2

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of inner packaging bag
Piece(s) / inner box 100
Piece(s) / case 1000
Piece(s) / pallet 28000
Case volume 0,0404 cbm
Weight of product 0,005 kg
Weight of case 5,6 kg
EAN of inner box 4038917153645
EAN case 4038917038225

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT