



### Product description

Order number	71.9923.212
Product description	Weigh tray, 70 ml, (LxW): 72 x 72 mm, material: PVC, white, 100 piece(s)/bag

### Size

Sample volume	70 ml
Max. volume	100 ml
Width of product	72 mm
Length of product	72 mm

### Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	white

### Purity & certification

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	70000
Depth of case	428 mm
Width of case	234 mm
Height of case	358 mm
Case volume	0.0359 cbm
Weight of product	0.0024 kg
Weight of case	5.6 kg
EAN of inner box	4038917153669
EAN case	4038917038249

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

