



Product description

Order number 71.9923.212
Product description Weigh tray, 70 ml, (LxW): 72 x 72 mm, material: PVC, white, 100 piece(s)/bag

Size

Sample volume 70 ml
Max. volume 100 ml
Width of product 72 mm
Length of product 72 mm

Material & colours

Product material Polyvinylchloride (PVC)
Colour of product white

Purity & certification

Product category not regulated
CE certified no CE mark
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	70000
Depth of case	355 mm
Width of case	425 mm
Height of case	235 mm
Case volume	0.0355 cbm
Weight of product	0.0024 kg
Weight of case	5.6 kg
EAN of inner box	4038917153669
EAN case	4038917038249

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

