



Product description

| | |
|---------------------|-----------------------------------------------------------------------|
| Order number | 81.970 |
| Product description | Stirring rod, length: 120 mm, material: PP, white, 1,000 piece(s)/bag |

Size

| | |
|-------------------|--------|
| Diameter | 3 mm |
| Length of product | 120 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | white |

Purity & certification

| | |
|-------------------|-------------------------|
| Product category | no medical device IVD |
| Latex-free | yes |
| Safety data sheet | no |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 10000 |
| Piece(s) / pallet | 400000 |
| Depth of case | 438 mm |
| Width of case | 268 mm |
| Height of case | 155 mm |
| Case volume | 0.0182 cbm |
| Weight of product | 0.0006 kg |
| Weight of case | 6.23 kg |
| EAN of inner box | 4038917269186 |
| EAN case | 4038917003506 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

