



Product description

Order number	84.1027.010
Product description	Adhesive label, (LxW): 62 x 40 mm, material: paper, white, label layout: Name + date, print: white/black, 5,000 piece(s)/roll

Product characteristics

Colour of print/label	white/black
-----------------------	-------------

Size

Width of product	40 mm
Length of product	62 mm

Material & colours

Product material	paper
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	roll
Piece(s) / inner box	5000
Piece (s) / outer case	20000
Piece(s) / pallet	1280000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0002 kg
Weight of case	4.44 kg
EAN of inner box	4038917124706
EAN case	4038917019354

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

