



Product description

Order number	84.2620.001
Product description	Adhesive label, (LxW): 33 x 12 mm, material: paper, white, 4,000 piece(s)/roll

Size

Width of product	12 mm
Length of product	33 mm

Material & colours

Product material	paper
Colour of product	white

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	4000
Type of smallest subpackaging	roll
Piece(s) / inner box	4000
Piece (s) / outer case	8000
Piece(s) / pallet	2592000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0001 kg
Weight of case	1.16 kg
EAN of inner box	4038917416948
EAN case	4038917416931

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

