



Product description

Order number 84.2620.002
Product description Adhesive label, (LxW): 42 x 12 mm, material: paper, white, 4,000 piece(s)/roll

Size

Width of product 12 mm
Length of product 42 mm

Material & colours

Product material paper
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 4000
Type of smallest subpackaging roll
Piece(s) / inner box 4000
Piece (s) / outer case 8000
Piece(s) / pallet 2592000
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0.0051 cbm
Weight of product 0.0001 kg
Weight of case 1.2 kg
EAN of inner box 4038917416917
EAN case 4038917416900

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

