



Product Specification
Drop collector, length: 112 mm, PP, tube 75 x 13 mm with funnel cup

Page 1



Product description

Order number	86.290
Product description	Drop collector, length: 112 mm, material: PP, tube 75 x 13 mm with funnel cup, funnel: HD-PE, for extracting drops for salmonella diagnosis, 250 piece(s)/bag

Product characteristics

Design	tube 75 x 13 mm with funnel cup
Application	for extracting drops for salmonella diagnosis

Size

Diameter of opening	40 mm
Length of product	112 mm

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0095 kg
Weight of case	5.45 kg
EAN of inner box	4038917081085
EAN case	4038917035620

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

