



Product Specification
Drop collector, length: 112 mm, PP, tube with funnel cup



Product description

Order number	86.290
Product description	Drop collector, length: 112 mm, material: PP, tube with funnel cup, funnel: HD-PE, for extracting drops for salmonella diagnosis, 250 piece(s)/bag

Product characteristics

Design	tube with funnel cup
Application	for extracting drops for salmonella diagnosis

Size

Diameter of opening	40 mm
Length of product	112 mm

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Drop collector, length: 112 mm, PP, tube with funnel cup

Page 2

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Safety data sheet	nein
Batched	no

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0,0657 cbm
Weight of product	0,0095 kg
Weight of case	5,45 kg
EAN of inner box	4038917081085
EAN case	4038917035620

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-11-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT