



### **Product description**

Order number	86.9924.531
Product description	Stomacher bag, (LxW): 300 x 170 mm, material: LD-PE, transparent, Food testing laboratories

### **Product characteristics**

Application	Food testing laboratories
Example of use	Prolonged massage in a Stomacher® machine

### **Size**

Width of product	170 mm
Length of product	300 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



### **Purity & certification**

Latex-free	yes
Pyrogen-free	no
Safety data sheet	nein

### **Packaging**

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	80000
Depth of case	458 mm
Width of case	238 mm
Height of case	373 mm
Case volume	0,0408 cbm
Weight of case	9,3 kg
EAN of inner box	4038917155380
EAN case	4038917398442

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

