



### **Product description**

Order number 91.2620.020  
Product description Thermo transfer foil, black, 1 piece(s)/roll

### **Material & colours**

Colour of product black

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 2  
Type of smallest subpackaging roll  
Piece(s) / inner box 1  
Piece (s) / outer case 2  
Piece(s) / pallet 648  
Depth of case 260 mm  
Width of case 183 mm  
Height of case 85 mm  
Case volume 0.004 cbm  
Weight of product 0.209 kg  
EAN of inner box 4038917899703  
EAN case 4038917899741

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

