



### **Product description**

Order number	95.813.001
Product description	Dispenser box, material: PS, lid assembled

### **Product characteristics**

Base shape	flat base
------------	-----------

### **Size**

Diameter of base	81 mm
Diameter of lid	83 mm
Diameter of opening	81 mm
Length of product	110 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	white
Closure material	Polystyrene (PS)
Colour of cap	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	35
Piece(s) / inner box	35
Piece (s) / outer case	35
Piece(s) / pallet	1260
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0623 kg
Weight of case	2.77 kg
EAN of inner box	4038917121644
EAN case	4038917074797

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

