



### **Product description**

Order number	95.950
Product description	Felt marker, black, alcohol-resistant, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

### **Product characteristics**

Label/ Print	with print
Colour of print/label	black
Example of use	for labelling plastic tubes and reaction tubes

### **Size**

Diameter	13 mm
Length of product	133 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	white/black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      10  
Type of inner packaging                  case  
Piece(s) / inner box                      10  
Piece(s) / case                              500  
Piece(s) / pallet                          45000  
Case volume                                0,0176 cbm  
Weight of product                        0,0161 kg  
Weight of case                              9,1 kg  
EAN of inner box                         4038917380355  
EAN case                                      4038917380348

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

