



Product description

Order number	95.950
Product description	Felt marker, black, alcohol-resistant, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

Product characteristics

Label/ Print	with print
Colour of print/label	black
Example of use	for labelling plastic tubes and reaction tubes

Size

Diameter	13 mm
Length of product	133 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	white/black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 10
Type of inner packaging case
Piece(s) / inner box 10
Piece(s) / case 500
Piece(s) / pallet 45000
Case volume 0,0176 cbm
Weight of product 0,0161 kg
Weight of case 9,1 kg
EAN of inner box 4038917380355
EAN case 4038917380348

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

