



Product description

Order number	95.955
Product description	Felt marker, green, waterproof, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

Product characteristics

Colour of print/label	white
Example of use	for labelling plastic tubes and reaction tubes

Size

Diameter	10 mm
Length of product	142 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Felt marker, green, waterproof

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 10
Type of smallest case
subpackaging
Piece(s) / inner box 10
Piece (s) / outer case 1000
Piece(s) / pallet 48000
Depth of case 300 mm
Width of case 291 mm
Height of case 392 mm
Case volume 0.0342 cbm
Weight of product 0.009 kg
Weight of case 10 kg
EAN of inner box 4038917084642
EAN case 4038917003414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

