



Product description

Order number 95.1760.100
Product description Rubber foot, grey

Product characteristics

Suitable for pipette tip boxes

Material & colours

Product material Thermoplastic elastomers (TPE)
Colour of product grey

Purity & certification

Product in accordance with MDD or IVDD no medical device | IVD
Lot or serial number no

Packaging

Minimum order qty. 25
Piece(s) / inner box 25
Piece(s) / case 2000
Piece(s) / pallet 648000
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0,0051 cbm
Weight of product 0,0001 kg
Weight of case 0,45 kg
EAN of inner box 4038917522670
EAN case 4038917522663

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

