



Product description

Order number	95.1763
Product description	Empty packaging, grey

Size

Width of product	128 mm
Height of product	112 mm
Length of product	98 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	grey
Colour of cap	transparent

Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	8
Piece(s) / inner box	8
Piece(s) / case	8
Piece(s) / pallet	624
Case volume	0,0148 cbm
Weight of product	0,1176 kg
Weight of case	1,27 kg
EAN of inner box	4038917087612
EAN case	4038917080378

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

