



Product description

Order number	70.1130.100
Product description	Pipette tip, working volume: 10 µl, transparent, filling level rings, suitable for Eppendorf, Gilson, Finnpipette, Brand and Socorex and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpipette, Brand and Socorex and products of identical design
Volume of work	10 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

Size

Max. volume	10 µl
Length of product	31 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1920
Type of inner packaging box
Piece(s) / inner box 96
Piece(s) / case 1920
Piece(s) / pallet 138240
Depth of case 358 mm
Width of case 323 mm
Height of case 140 mm
Case volume 0.0162 cbm
Weight of product 0.0011 kg
Weight of case 2.33 kg
EAN of inner box 4038917125864
EAN case 4038917020312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

