



Product description

Order number	70.1130.100
Product description	Pipette tip, working volume: 10 µl, transparent, filling level rings, suitable for Eppendorf, Gilson, Finnpipette, Brand and Socorex and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpipette, Brand and Socorex and products of identical design
Volume of work	10 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	10 µl
Length of product	31 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Packaging

Minimum order qty.	1920
Type of inner packaging	box
Piece(s) / inner box	96
Piece(s) / case	1920
Piece(s) / pallet	138240
Depth of case	358 mm
Width of case	323 mm
Height of case	140 mm
Case volume	0,0162 cbm
Weight of product	0,0011 kg
Weight of case	2,33 kg
EAN of inner box	4038917125864
EAN case	4038917020312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

