



### **Product description**

Order number	70.1183.001
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Gilson, 250 piece(s)/bag

### **Product characteristics**

Graduation	yes
Suitable for	Gilson
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	yes

### **Packaging**

Minimum order qty.	2500
Type of inner packaging	bag
Piece(s) / inner box	250
Piece(s) / case	2500
Piece(s) / pallet	25000
Case volume	0,0737 cbm
Weight of product	0,0026 kg
Weight of case	7,3 kg
EAN of inner box	4038917126441
EAN case	4038917080446

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

