



### Product description

Order number	70.1183.002
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Eppendorf, 250 piece(s)/bag

### Product characteristics

Graduation	yes
Suitable for	Eppendorf
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

### Size

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material                      Polypropylene (PP)  
Colour of product                      transparent

### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      2500  
Type of inner packaging                      bag  
Piece(s) / inner box                      250  
Piece(s) / case                              2500  
Piece(s) / pallet                              25000  
Case volume                                0,0966 cbm  
Weight of product                        0,0026 kg  
Weight of case                              7,3 kg  
EAN of inner box                        4038917126458  
EAN case                                      4038917080453

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

