



Product description

| | |
|---------------------|---|
| Order number | 70.1183.101 |
| Product description | Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Gilson, 50 piece(s)/box |

Product characteristics

| | |
|------------------------|--------|
| Graduation | yes |
| Suitable for | Gilson |
| Volume of work | 5 ml |
| Filter | no |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|--------|
| Max. volume | 5 ml |
| Length of product | 123 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | box |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 16500 |
| Depth of case | 165 mm |
| Width of case | 311 mm |
| Height of case | 610 mm |
| Case volume | 0.0313 cbm |
| Weight of product | 0.0065 kg |
| Weight of case | 3.7 kg |
| EAN of inner box | 4038917126465 |
| EAN case | 4038917080460 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

